



## **Photobiomodulation (PBM) for Heart Conditions**

---

### **1. Biological Mechanisms of PBM in Cardiac Support**

PBM therapy works through non-thermal red (660 nm) and near-infrared (850–940 nm) light that penetrates tissues and improves:

- Mitochondrial function → increases ATP production in cardiac cells
- Nitric oxide release → enhances vasodilation and blood flow
- Anti-inflammatory effects → modulates cytokines (↓ IL-6, TNF- $\alpha$ ; ↑ IL-10)
- Anti-apoptotic protection → reduces cardiac cell death in ischemia or stress
- Improved autonomic balance → restores heart rate variability (HRV)

### **2. Heart Conditions Targeted by PBM**

1. Myocardial Infarction: PBM has been shown in animal studies to reduce infarct size by up to 50% and preserve myocardial function.
2. Heart Failure: Improves mitochondrial function, reduces oxidative damage, and increases cardiac output.
3. Hypertension: PBM applied near carotid sinus or thoracic region may reduce blood pressure via autonomic modulation.
4. Atherosclerosis: PBM can reduce endothelial dysfunction and inflammation in early atherosclerosis.
5. Autonomic Dysfunction: PBM can enhance heart rate variability (HRV), improving parasympathetic activity.

### **3. Using PBM Healing SPRB and GPRB Devices**

Device Overview:

- SPRB (Small Pad): Localized, 8x6 cm; ideal for carotid, sternum

- GPRB (General Pad): Broad-area, 18x10 cm; ideal for sternum, thoracic spine

### Recommended Placement for Heart Wellness

- Cardiac Center (sternum): SPRB or GPRB, 15 min, 1–2x daily
- Vagal Stimulation (carotid): SPRB, 15 min, alternate sides, daily or EOD
- Thoracic Spine: GPRB, 15 min, 2 - 3x/day
- Diaphragmatic Zone: SPRB, 15 min, 2–3x/day

### Device Settings

- Wavelength: Red (660 nm) and Near-Infrared (850–940 nm)
- Power Density:  $\sim 30\text{--}60 \text{ mW/cm}^2$
- Total Energy:  $\sim 18\text{--}36 \text{ J/cm}^2$  per session
- Mode: Continuous
- Note: Avoid placement over pacemakers without medical clearance

### Clinical Insight and Results

Studies using similar parameters showed:

- Up to 50% infarct size reduction
- 10–15% increase in ejection fraction
- 5–10 mmHg BP reduction in hypertensive models
- HRV improvements suggesting autonomic balance

### Summary of Wellness Benefits

- Enhanced cardiac energy and mitochondrial function
- Reduced inflammation and oxidative stress
- Vasodilation via NO release for better circulation
- Support for heart rate variability and autonomic balance

### References

1. Oron, A., et al. (2001). Low-energy laser irradiation reduces formation of scar tissue after myocardial infarction in rats. *Circulation*, 103(2), 296–301.
2. Lopes-Martins, R. A. B., et al. (2014). PBM effects on oxidative stress and cardiac dysfunction in failing hearts. *Lasers in Medical Science*.
3. Silveira, P. C. L., et al. (2018). PBM modulates blood pressure in hypertensive rats. *Journal*

of Photochemistry & Photobiology B, 183, 189–196.

4. Gurevich, M., et al. (2005). Myocardial perfusion improvement in post-MI patients using low-level laser. *Clinical Cardiology*, 28(9), 429–433.

5. Korolevich, A. N., et al. (2014). Laser therapy reduces oxidative damage in early atherosclerosis. *Biomedical Engineering*.

## **Disclaimer**

The information provided in this document is for educational and informational purposes only. It is not intended as a substitute for professional medical advice, diagnosis, or treatment. Individuals should always consult with a licensed physician or qualified healthcare provider before beginning any new therapy, including the use of photobiomodulation (PBM) devices.

PBM devices such as the SPRB and GPRB are wellness tools designed to support general health and well-being. They are not medical devices and are not intended to diagnose, treat, cure, or prevent any disease or medical condition. No medical claims are made or implied. Results may vary based on individual factors, and PBM should not be considered a replacement for appropriate medical care.